**U L T R A S O N I C   I N S E R T S**  
**U L T R A S C H A L L   E I N S Ä T Z E**

---

Piezo Tips/ <i>Ultraschallspitzen</i> .....	C1
Design Options/ <i>Design Optionen</i> .....	C2
General Debridement .....	C3
Periodontal Debridement.....	C6
Reference/ <i>Index</i> .....	C7

# Just what I expect from Hu-Friedy.

## New **Symmetry** S-Series Piezo Scaling Tips

Hu-Friedy – the world leader in scalers and curettes – introduces Symmetry Piezo Tips. Our special blend of stainless steel works in harmony with your generator's sonic energy to optimize vibration for superior scaling efficiency and unmatched durability. S-Series tips are compatible with all Satelec generators.

*Generating Perfect Harmony.*

## the **Guardian** Tip Carrier

Symmetry S-Series Piezo Tips come with the Guardian Tip Carrier – an innovative torque limiting wrench/tip carrier designed for remarkably safe and easy connection to the handpiece. The exclusive power range indicator alerts clinicians to required power setting for optimal tip performance.

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x\* See Section A for exact product code number for cassette color desired.

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The GUARDIAN™ Tip Carrier Symmetry S-Series Piezo Tips come with the Guardian Tip Carrier – an innovative torque limiting wrench/tip carrier designed for remarkably safe and easy connection to the handpiece. The exclusive power range indicator alerts clinicians to required power setting for optimal tip performance.

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**Furcation PLUS™ Designs**

**Efficient at Lower Power Ranges**

Allows you to turn down the power to increase patient comfort and still achieve efficient, thorough scaling.

**Color-Coded Handles**

Colorful resin and silicone handles help instantly identify tip design and increase tactile sensitivity.

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25kHz Part Code  
30kHz Part Code

**SATIN SWIVEL™ DESIGNS Revolutionary SWIVEL™ feature for smooth scaling from arch to arch or tooth to tooth**

Experience the maneuverability and adaptability you want in power scaling with the exclusive SATIN SWIVEL™ inserts. With a slight rolling motion of your fingertips, you can eliminate the old "start, stop, remove and readjust" routine required with traditional inserts. SATIN SWIVEL™ inserts also reduce the torque from the handpiece cord, increasing overall balance for hassle-free scaling.

The grip – up to 26% larger than most standard inserts – reduces finger pinching, resulting in less hand fatigue\* and enhanced tactile sensitivity. The lightly textured silicone grips are color-coded for easy tip identification and provide a secure, comfortable grasp.

\*Source: Gerwatowski, L.J., McFall, D.B., Stach, D.: Carpal Tunnel Syndrome; Risk Factors and Preventive Strategies for the Dental Hygienist. Journal of Dental Hygiene, February 1992.

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\*Quelle: Gerwatowski, L.J., McFall, D.B., Stach, D.: Carpal Tunnel Syndrome; Risk Factors and Preventive Strategies for the Dental Hygienist. Journal of Dental Hygiene, February 1992.

**SWIVEL DIRECT FLOW ULTRASONIC INSERTS/SWIVEL DIRECT FLOW ULTRASCHALL EINSÄTZE**

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\*\*] shown in = 30 kHz size  
\*\*] Abbildung = 30 kHz

**DID YOU KNOW?**

Hu-Friedy inserts are compatible with all Cavitron style magnetostrictive scaling units.

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GENERAL DEBRIDEMENT

CE0120

STREAMLINE® DESIGNS / STREAMLINE® DESIGNS

#10 Universal

UI1025K

UI1030K

25kHz Part Code  
30kHz Part Code

#100 Thin/dünn

UI25K100S

UI30K100S

#1000 Triple Bend/3fach

UI100025K

UI100030K

#3 Beavertail

UI325K

UI330K

ORIGINAL PROPHY DESIGNS / ORIGINAL PROPHY DESIGNS

Original Propphy Style

UI25KP10

UI30KP10

Original Propphy PLUS™

UI25KP10P

UI30KP10P

ORIGINAL PROPHY DESIGNS

The standard diameter Propphy style tip is the original ultrasonic insert design, with a finely beveled external water delivery tube. The PLUS™ design features a fixed water delivery tube in a lightweight resin handle.

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STREAMLINE® DESIGNS

Efficient at low power settings, internal flow water delivery

Routine scaling is more efficient with our Streamline® inserts, even at low power settings. Designed to remove general deposits in less time, the water is delivered directly from the tip base, eliminating the need for an external water tube. Resin handles are color-coded to help instantly identify tip style and increase tactile sensitivity. The #3 Beavertail design is excellent for removal of gross calculus and orthodontic cement.

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DID YOU KNOW?

The SATIN SWIVEL™ and the Streamline® #100 Thin inserts, which feature a 40% thinner tip than the standard #10 insert, are ideal for light deposits. Recommended for use on LOW TO MEDIUM power only.

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**PERIODONTAL DEBRIDEMENT**

CE 0120

**AFTER FIVE® PLUS™ DESIGNS/AFTER FIVE® PLUS™ DESIGNS**

Left/links

25kHz Part Code UI25KSF10L

30kHz Part Code

UI30KSF10L

Straight/gerade

25kHz Part Code UI25KSF10S

30kHz Part Code UI30KSF10S

Right/rechts

25kHz Part Code UI25KSF10R

30kHz Part Code UI30KSF10R

**AFTER FIVE® PLUS™ DESIGNS**

**Fixed water delivery; thinner tip design**

Reach deep into pockets and achieve efficient subgingival scaling with these thin tips, exclusively from Hu-Friedy. Probe-size tips (same tip design as the After Five® inserts) offer smooth, deep access into periodontal pockets. A fixed water tube provides hassle-free water delivery. Patient comfort is maximized since designs are efficient at low power settings. Our resin handle is lightweight and color-coded for instant tip identification, increased tactile sensitivity, and reduced operator fatigue.

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**AFTER FIVE® DESIGNS/AFTER FIVE® DESIGNS**

Left/links

25kHz Part Code UI25KSL10L

30kHz Part Code UI30KSL10L

Straight/gerade

25kHz Part Code UI25KSL10S

30kHz Part Code UI30KSL10S

Right/rechts

25kHz Part Code UI25KSL10R

30kHz Part Code UI30KSL10R

PERIODONTAL DEBRIDEMENT

CE0120

FURCATION PLUS™ DESIGNS / FURCATION PLUS™ DESIGNS

Furcation PLUS™ Left/links

UI25KFPL  
UI30KFPL

Furcation PLUS™ Straight/gerade

UI25KFPS  
UI30KFPS

Furcation PLUS™ Right/rechts

UI25KFPR  
UI30KFPR

0.8mm ball end  
Rundes 0,8 mm-Ende

FURCATION DESIGNS

0.8mm ball end design adapts to furcations

Ball end offers greater tip surface area for more thorough debridement in this difficult to reach area. Also ideal for developmental concavities. The slim design allows greater access with less tissue distension and the water delivery tube is positioned for focused lavage. In the PLUS™ design, the water tube is fixed within a lightweight resin handle for hassle-free water delivery. These unique Hu-Friedy Furcation designs are available in right, left, and straight patterns for a full range of access.

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**CLEANING**

To clean, rinse thoroughly or, if you choose, fully immerse in a mild ultrasonic cleaning solution (pH 7-10.5) like Hu-Friedy's IMS® Daily Clean or Enzymax®. Place in ultrasonic cleaning unit for 7-10 minutes or 16-20 minutes if the insert is in a cassette. Thoroughly rinse the insert with tap water and dry completely prior to sterilization.

**STERILIZATION**

Sterilize Hu-Friedy Ultrasonic Inserts in an instrument cassette (i.e. IMS instrument cassette) or paper pouch in a steam autoclave.

Chemical vapor sterilization is unacceptable for SATIN SWIVEL™ inserts. Rapid heat transfer sterilization and cold liquid/chemical sterilants are not recommended and may void warranty. Avoid contact with corrosive chemicals such as chlorides, sulfates and surface disinfectants. Do not expose to solutions containing phenols.

**CARE AND MAINTENANCE**

To increase the o-ring life of Hu-Friedy Ultrasonic Inserts and to maintain a smooth swivel function, take the following steps before placing the insert into the handpiece. Fill entire handpiece with water. Lubricate the o-ring with water. Gently twist the insert into the handpiece until completely seated. For Satin Swivel™ inserts, activate the foot switch and rotate insert 2-3 times.

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
Scaling efficiency can significantly diminish with worn insert tips. Inferior performance and poor water delivery can result from worn, damaged, bent, or altered tips.

For optimal scaling efficiency, we recommend checking inserts monthly for signs of wear.

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**DID YOU KNOW?**


All thin inserts (SATIN SWIVEL™ #10, Streamline® #100, After Five®, After Five® PLUS™, Furcation, Furcation PLUS™) should be used on LOW TO MEDIUM power settings only.

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LIGHT DEBRIS REMOVAL / DEPLAQUING

- thin, perio design inserts  
SATIN SWIVEL™ #100, Streamline® #100,  
After Five®, After Five® PLUS™
- low to medium power setting only

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
SITE-SPECIFIC DEBRIS REMOVAL

- thin inserts with 0.8mm ball end  
Furcation, Furcation PLUS™
- low to medium power setting only

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DID YOU KNOW?

By varying the position of the insert tip on the tooth you can increase deposit removal efficiency (concave and convex surfaces) or increase patient comfort (lateral surfaces).

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REFERENCE SECTION/INDEX

CROSS REFERENCE GUIDE/INDEX

HU-FRIEDY ULTRASONIC INSERTS HU-FRIEDY ULTRASCHALL-EINSÄTZE				DENTSPLY CAVITRON* INSERTS DENTSPLY CAVITRON* EINSÄTZE		
Tip Style Spitze	Handle Griff	25kHz Part Code	30kHz Part Code	Tip Style Spitze	25kHz Handle Griff	30kHz Handle Griff
<b>SWIVEL™ DIRECT FLOW</b>						
#10 Universal	lavender	UI25SD10	UI30SD10			
#100 Thin	black	UI25SD100	UI30SD100			
#After Five™ Left	teal	UI25SDL	UI30SDL			
#After Five™ Straight	darkblue	UI25SDS	UI30SDS			
#After Five™ Right	red	UI25SDR	UI30SDR			
<b>SATIN SWIVEL™</b>						
#10 Universal	lavender	UI25SS10	UI30SS10			
#100 Thin	black	UI25SS100	UI30SS100			
#1000 Triple Bend	orange	UI25SS1000	UI30SS1000			
<b>Streamline®</b>				<b>TFI* Thru Flow*</b>		
#10 Universal	lavender	UI1025K	UI1030K	TFI #10 Straight	grey	black
#100 Thin	black	UI25K100S	UI30K100S			
#1000 Triple Bend	orange	UI100025K	UI100030K	TFI #1000 Straight Triple	grey	black
#3 Beavertail	yellow	UI325K	UI330K	Bend TFI #3 Thin Beavertail	grey	black
<b>Original Prophy Designs</b>				<b>Classic "P"*</b>		
#10 Universal	metal	UI25KP10	UI30KP10	25K P-10 Classic "P"	metal	metal
#10 Universal PLUS™	green	UI25KP10P	UI30KP10P			
<b>After Five® – Periodontal Debridement</b>				<b>SLI* Slimline*</b>		
After Five® Straight	metal	UI25KSL10S	UI30KSL10S	SLI #10 Straight/P-12	metal	metal
After Five® Right	metal	UI25KSL10R	UI30KSL10R	SLI #10 Right/EWP-12R	metal	metal
After Five® Left	metal	UI25KSL10L	UI30KSL10L	SLI #10 Left/EWP12L	metal	metal
				SLI Assorted 3-pack	metal	metal
<b>After Five® PLUS™ – Periodontal Debridement</b>				<b>FSI* Focused Spray* and FSI-SLI Focused Spray Slimline</b>		
After Five® PLUS™ Straight	dark blue	UI25KSF10S	UI30KSF10S	Refer to Hu-Friedy's After Five® PLUS™ inserts for an appropriate application. <i>Hinweise zur fachgerechten Anwendung finden Sie im Abschnitt mit den After Five® PLUS™ Einsätzen von Hu-Friedy.</i>		
After Five® PLUS™ Right	red	UI25KSF10R	UI30KSF10R			
After Five® PLUS™ Left	teal	UI25KSF10L	UI30KSF10L			
<b>Furcation and Furcation PLUS™ – Periodontal Debridement</b>						
Furcation Straight	metal	UI25KFS	UI30KFS			
Furcation Right	metal	UI25KFR	UI30KFR			
Furcation Left	metal	UI25KFL	UI30KFL			
Furcation PLUS™ Straight	cream	UI25KFPS	UI30KFPS			
Furcation PLUS™ Right	brown	UI25KFPR	UI30KFPR			
Furcation PLUS™ Left	light blue	UI25KFPL	UI30KFPL			

\*CAVITRON is a trademark of DENTSPLY International, Inc. and is registered with the U.S. Patent and Trademark Office.  
SLIMLINE, THRU FLOW, CLASSIC "P", TFI, SLI, FSI, FOCUSED SPRAY are trademarks of DENTSPLY International, Inc.  
\*CAVITRON ist ein Warenzeichen von DENTSPLY International, Inc. und beim USPTO („U.S. Patent and Trademark Office“) eingetragen.  
SLIMLINE, THRU FLOW, CLASSIC "P", TFI, SLI, FSI, FOCUSED SPRAY sind Warenzeichen von DENTSPLY International, Inc.